EXAMPLE FORM EQUIPMENT PERFORMANCE EVALUATION DENTAL RADIOGRAPHIC UNIT 25 TAC §289.232(i)(7)(A) Facility Name:					
	ce Company:Service Company Registration No.:				
Survey Instrument Used:	ey Instrument Used:			nsor/detector	
Instrument Calibration/ Intercomparison	ı Date:	Technician Signa	ture:		
X-RAY UN	NIT IDENTIFICATION (F	FROM CONTROL F	PANEL)		
Manufacturer:			Location/Room:		
Model No.:	Serial N	No.:			
	TIMER ACCU	RACY	(	) Pass ( ) Fail	
Regulation –The accuracy of the timer sha obtainable, the timer accuracy shall be $\pm 1$ §289.232(i)(6)(H)( <i>i</i> ) )					
SELECT ONE: D Manufacturer's spec	ifications which are		(	)R	
$\Box$ Intraoral units ±10 %	% tolerance with testing perf	formed at 0.5 second (	500 milliseconds)		
$\Box$ Panoramic units ±10	0% tolerance with testing per	rformed at	seconds		
Document time used for testing:	msec; OR	pulses;	OR	seconds	
Perform four measurements: 1	2	3	4		
	EVDOCUDE DEDDOD		(	) Desg. ( ) Feil	
Regulation – When all technique factors an systems shall not exceed 0.05. (See 25 Ta		ent of variation of exp		) Pass ( ) Fail	
Technique factors selected:mA_	kVptim				
Perform four measurements:			ndard deviation of t		
1mR	$-2 \frac{1}{2}$	$X_i$ = ith observati	on in sample		
<sup>2.</sup> mR $C = \frac{s}{m} = \frac{l}{m}$	$\sum_{i=1}^{n} \frac{(X_i - \overline{X})^2}{n - 1} \Big ^{1/2}$	N = number of ol	oservations in samp	le	
3. <u>m</u> R X X 4mR	$\sum_{i=1}^{n-1}$	Coefficient of var (Must not exceed .			

KVP TEST( ) Pass( ) FailRegulation – If the registrant possesses the manufacturer's kilovolt peak specifications, the radiation machine shall meetthose specifications. Otherwise, the measured kVp shall be accurate to within $\pm 10$ percent of the indicated setting at no less thanthree points over the usual operating range of the machine. (For units with fewer than three fixed kVp settings, the units shall bechecked at those settings.) (See 25 TAC §289.232(i)(6)(J)						
SELECT BELOW:		((Measured kVp - Ind	dicated kVp) ÷Inc	dicated kVp) x $100 =$	% of Deviation	
☐ Manufacturer's	specifications which	are		OR		
$\Box \pm 10\%$ of the indi	cated setting used					
		1. Indicated kVp 2. Indicated kVp			iation iation	% %
		3. Indicated kVp 4. Indicated kVp	Measured k	VpDevi	ationation	/0 % %
		TUBE ST	ABILITY		( ) Pass	( ) Fail
Regulation –The tube shall remain physically stable during exposures. In cases where tubes are designed to move during exposure the registrant shall assure proper and free movement of the unit. (See 25 TAC § 289.232(i)(6)(K))						
For <b>intraoral machine</b> Tube stable in all orien		es ( ) No ( )	-	and cephalometric mad where designed? Yes (	-	
		COLL	IMATION		( ) Pass	( ) Fail
Regulation-25 TAC §289.232 (i)(6)(L) Field limitation shall meet the requirements of 25 TAC §289.232(i)(11)(B) and 25 TAC §289.232(i)(12)						
Intraoral:		o skin distance (SSD) if the minimum SSD is				
Panoramic:	X-ray field size:	size: (Circle inch or ce soment at image recepto	Trans	sversein/cm sversein/cm sversein/cm	Vertical	in/cm in/cm in/cm
	(Misalignment can	not exceed 0.0 inches in	n the transverse as	kis and 0.5 inches in t	he vertical axis.)	
Cephalometric, Digital Panoramic, CBCT	Source to image d Image Receptor, E Measured x-ray fie X-ray field misali	Detector, FOV, Diaphra eld size:		_in/cm _in/cm xin _in/cm xin _in/cm xin	n/cm	
Alternative Results (please see attachments for documentation)						
		ENTRANCE	E EXPOSURE (E	 E)	()Pass	( ) Fail
Regulation – The in-air exposure limits for an average adult intraoral bite wing examination shall not exceed <b>450 mR for 60 kVp and above; 600 mR for less than 60 kVp.</b> (See 25 TAC  289.232(i)(6)(M) ) (See page 5 for instructions.)						
Technique factors, for intraoral bite wing examination only, selected: kVp mA(s) time   Source to Skin Distance (SSD)in/cm Source to Detector Distance (SDD)in/cm						
Tip of cone is positioned ½ inch or less from surface of instrument housing or detector. Yes ( ) No( )						
Calculated Meas Detector Meas.	surement <b>OR</b> D mR	irect Measurement		EE:	mR	

## EXPOSURE REPRODUCIBILITY CALCULATIONS **EXAMPLE**

$$C = \frac{s}{\overline{X}} = \frac{l}{\overline{X}} \left[ \sum_{i=1}^{n} \frac{(X_i - \overline{X})^2}{n - l} \right]^{1/2}$$

Where:

C = coefficient of variations = estimated standard deviation of the population  $\overline{X}$  = mean value of observations in sample  $X_i = ith \ observatio \ n \ in \ sample$ n = number of observations in sample

In this example, the exposures are considered to be reproducible.

Example:

The four (n) exposures ( $\chi_i$ ) measured 409 mR, 387 mR, 391 mR, and 410 mR.

Determine the mean value  $(\overline{X})$  of the four exposures taken. STEP 1

 $(409 mR + 387 mR + 391 mR + 410 mR) \div 4 = 399.25 mR$ 

STEP 2	Find the difference between each exposure and the mean value (disregard sign).			
	409.00 mR	387.00 mR	391.00 mR	410.00 mR
	<u>-399.25 mR</u>	<u>-399.25 mR</u>	<u>-399.25 mR</u>	<u>-399.25 mR</u>
	9.75 mR	12.25 mR	8.25 mR	10.75 mR

STEP 3 Square each of the differences

> $12.25^2 = 150.06$  $9.75^2 = 95.06$  $10.75^2 = 115.56$  $8.23^2 = 67.73$

STEP 4	Divide each number by 3 $(n-1)$ and add the results
	95.06 - 3 = 31.69 150.06 - 3 = 50.02 68.06 - 3 = 22.69 $115.56 - 3 = \underline{38.52}$ 142.92
STEP 5	For s, determine the square root of the above number
	$\sqrt{142.92} = 11.96$
STEP 6	Divide s by the mean value $(\bar{x})$
	11.96 - 399.25 = .0299 = C = the coefficient of variation
STEP 7	If C=0.05 or less, the exposures are considered to be reproducible

## ENTRANCE EXPOSURE (EE) CALCULATIONS EXAMPLE – (DENTAL – BITE WING)

## SENSOR/DECTECTOR ENCLOSED WITHIN A HOUSING

Where:

EE = entrance exposure mR (MEASURED) = indicated exposure SDD = source (target) to enclosed detector distanceSSD = source to surface (skin) distance

## For this example only: SDD = 24 cm. SSD = 21 cm.

Detector = 2cm. below housing surface Measurement = 374 mR (Verify that all units of measurements are the same – either inches or centimeters)

STEP 1	<b>Set the control panel</b> with the technique factors (kVp, mA, time) that the facility uses for an average person intraoral bite wing examination.
STEP 2	<b>Measure the SSD</b> . To obtain an accurate EE, have the operator position the tube over the detector housing as if they were positioning a patient. Measure from the target to the housing surface. (21 cm.)
STEP 3	<b>Measure the SDD</b> . Place the tip of the cone 1.5 cm. or less from the housing surface. Measure from the target to the surface of the housing and add the predetermined distance of the detector below the surface. $22 \text{ cm.}+2 \text{ cm.} = 24 \text{ cm.}$
STEP 4	Divide the SDD by the answer in Step 2. $24 \div 21 = 1.143$
STEP 5	Square the answer in Step 4. 1.143 x 1.143 = 1.307
STEP 6	Multiply the measurement reading in mR by the answer in Step 5. $374 \times 1.307 = 489 \text{ mR}$
STEP 7	Compare the answer in Step 6 to the regulatory limits of 450 mR when 60 kVp or above is used.